



An uplifting sound experience

Users recommend Interton Presto for sound quality and hearing in noise. Start enjoying clear and easy conversations.

Styles designed to delight you

With the small, comfortable designs and our latest technology, you'll look forward to putting on your Interton Presto hearing aids.

Greater ways to stay connected

Connect with Bluetooth® Low Energy (LE) Audio and Auracast broadcast audio, the new global Bluetooth experience and assistive listening technology.



Hear your best anywhere

Focus on the details you enjoy or who you want to listen to while hearing clearly all around you.

Enjoy hands-free calls

Stream hands-free calls from iPhone, iPad and compatible Android™ mobile devices* directly to your Interton Presto hearing aids.

* Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.
Bluetooth LE Audio works with Interton Presto and compatible mobile devices with Bluetooth 5.3 or later. See compatibility on interton.com/compatibility

What's your ideal style?

In the ear, behind the ear or custom made, there's a smaller-than-ever style that's great for you.

100% recyclable packaging

Your Interton Presto hearing aids come in packaging created from FSC certified paper and with 47% reduced CO2 emissions** as part of our focus on sustainability.

** Compared to previous hearing aid packaging from GN, CO2 reduction verified by Bureau Veritas.

Stylish wireless accessories

Designed for life's loud moments and made for Auracast™. Stream sound to boost your moments.

Entertainment and conversations

Hear TV clearly with the TV-Streamer+ and cut through the noise with the Multi-Mic+.

	Feature name	Presto 9	Presto 7	Presto 5	Presto 4
Audiological features	Number of channels	17	14	12	12
	M&RIE	✓	✓	✓	✓
	All access directionality ii	✓	✗	✗	✗
	Spatial sound	✓	✓	✓	✗
	Front focus	✓	✗	✗	✗
	Synchronised auto-steered directionality	✓	✓	✓	✓
	Environmental tuner ii	✓	✗	✗	✗
	Adaptive noise reduction	5 Settings (10 dB)	3 Settings (8 dB)	2 Settings (6 dB)	On / Off (3 dB)
	Microphone noise reduction	3 Settings	2 Settings	On / Off	On / Off
	Impulse noise reduction	3 Settings	3 Settings	On / Off	On / Off
	Adaptive wind noise reduction	3 Settings	2 Settings	On / Off	On / Off
	Frequency shifter	✓	✓	✓	✓
	Feedback manager plus (w/ music mode)	✓	✓	✓	✓
	Synchronised acclimatization manager	✓	✓	✓	✓
Functional features	Ear-to-ear communication	✓	✓	✓	✓
	Bluetooth® Auracast™	✓	✓	✓	✓
	Direct audio streaming (compatible ios** and android™ devices)	✓	✓	✓	✓
	Tv streamer 2, tv-streamer+, remote control, remote control 2, phone clip+	✓	✓	✓	✓
	Micro Mic, Multi Mic & Multi-Mic+	✓	✓	✓	✓
	Interton sound app	✓	✓	✓	✓
	Sound enhancer (interton sound app)	✓	✗	✗	✗
Online services	Interton remote fitting	✓	✓	✓	✓
	Interton remote fitting live	✓	✓	✓	✓
	Remote firmware updates	✓	✓	✓	✓
Fitting features	Fully flexible programs	4	4	4	4
	Pre-set feedback manager	✓	✓	✓	✓
	Datalogging	✓	✓	✓	✓
	Tinnitus sound generator	✓	✓	✓	✓
	Automatic receiver detection	✓	✓	✓	✓
	Wireless fitting with noahlink wireless	✓	✓	✓	✓



Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP. in ₹	
NEW		Presto Series 9			Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher					Chipset - C 6			
1	Presto 9 microRIE (2 Units + 1 Premium charger)	MM / LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	7,00,000 /-
2	Presto 9 microRIE (1 Unit + 1 Standard charger)	MM / LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	3,60,000 /-
3	Presto 9 HP BTE (2 Units + 1 Desktop charger)		Noahlink wireless	17	Li-ion	80	30 - 110	7,00,000 /-
4	Presto 9 HP BTE (1 Unit + 1 Desktop charger)		Noahlink wireless	17	Li-ion	80	30 - 110	3,58,000 /-
5	Presto 9 NP BTE (2 Units + 1 Desktop charger)		Noahlink wireless	17	Li-ion	67	20 - 100	7,00,000 /-
6	Presto 9 NP BTE (1 Unit + 1 Desktop charger)		Noahlink wireless	17	Li-ion	67	20 - 100	3,58,000 /-
7	Presto 9 ITC/ITE (2 Units + 1 Custom charger)	LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	50 / 60 / 70 / 78	20 - 110	7,00,000 /-
8	Presto 9 ITC/ITE (1 Unit + 1 Custom charger)	LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	50 / 60 / 70 / 78	20 - 110	3,60,000 /-
9	Presto 9 ITC	LP / MP / HP / UP	Noahlink wireless	17	312	51 / 60 / 69 / 77	20 - 110	3,50,000 /-
10	Presto 9 NWL CIC	LP / MP / HP	Flex cable / HI PRO	17	10A	51 / 62 / 69	20 - 105	3,50,000 /-

Presto Series 9 2 + 2 Years = 4 Years Total Warranty

NEW			Presto Series 7			Made for iPhone iPad iPod		Works with android		Made for Android™ ASHA		
Software - Interton fitting 2.0 & higher										Chipset - C6		
11	Presto 7 microRIE (2 Units + 1 Premium charger)	MM / LP / MP / HP / UP		Noahlink wireless	14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	4,00,000 /-			
12	Presto 7 microRIE (1 Unit + 1 Standard charger)	MM / LP / MP / HP / UP		Noahlink wireless	14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,10,000 /-			
13	Presto 7 HP BTE (2 Units + 1 Desktop charger)			Noahlink wireless	14	Li-ion	80	30 - 110	3,50,000 /-			
14	Presto 7 HP BTE (1 Unit + 1 Desktop charger)			Noahlink wireless	14	Li-ion	80	30 - 110	1,83,000 /-			
15	Presto 7 NP BTE (2 Units + 1 Desktop charger)			Noahlink wireless	14	Li-ion	67	20 - 100	3,50,000 /-			
16	Presto 7 NP BTE (1 Unit + 1 Desktop charger)			Noahlink wireless	14	Li-ion	67	20 - 100	1,83,000 /-			
17	Presto 7 ITC/ITE (2 Units + 1 Custom charger)	LP / MP / HP / UP		Noahlink wireless	14	Li-ion	50 / 60 / 70 / 78	20 - 110	3,50,000 /-			
18	Presto 7 ITC/ITE (1 Unit + 1 Custom charger)	LP / MP / HP / UP		Noahlink wireless	14	Li-ion	50 / 60 / 70 / 78	20 - 110	1,85,000 /-			
19	Presto 7 ITC	LP / MP / HP / UP		Noahlink wireless	14	312	51 / 60 / 69 / 77	20 - 110	1,75,000 /-			
20	Presto 7 NWL CIC	LP / MP / HP		Flex cable / HI PRO	14	10A	51 / 62 / 69	20 - 105	1,75,000 /-			

Presto Series 7 2 + 2 Years = 4 Years Total Warranty

Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP. in ₹	
NEW		Presto Series 5			Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher						Chipset - C 6		
21	Presto 5 microRIE (2 Units + 1 Premium charger)	MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,50,000 /-
22	Presto 5 microRIE (1 Unit + 1 Standard charger)	MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,35,000 /-
23	Presto 5 HP BTE (2 Units + 1 Desktop charger)		Noahlink wireless	12	Li-ion	80	30 - 110	2,00,000 /-
24	Presto 5 HP BTE (1 Unit + 1 Desktop charger)		Noahlink wireless	12	Li-ion	80	30 - 110	1,08,000 /-
25	Presto 5 NP BTE (2 Units + 1 Desktop charger)		Noahlink wireless	12	Li-ion	67	20 - 100	2,00,000 /-
26	Presto 5 NP BTE (1 Unit + 1 Desktop charger)		Noahlink wireless	12	Li-ion	67	20 - 100	1,08,000 /-
27	Presto 5 ITC/ITE (2 Units + 1 Custom charger)	LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	2,20,000 /-
28	Presto 5 ITC/ITE (1 Unit + 1 Custom charger)	LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	1,20,000 /-
29	Presto 5 ITC	LP / MP / HP / UP	Noahlink wireless	12	312	51 / 60 / 69 / 77	20 - 110	1,00,000 /-
30	Presto 5 NWL CIC	LP / MP / HP	Flex cable / HI PRO	12	10A	51 / 62 / 69	20 - 105	1,00,000 /-

Presto Series 5 2 + 1 Years = 3 Years Total Warranty

NEW		Presto Series 4			Made for iPhone iPad iPod		Works with android	Made for Android™ ASHA		
Software - Interton fitting 2.0 & higher								Chipset - C6		
31	Presto 4 microRIE (2 Units + 1 Premium charger)	MM / LP / MP / HP / UP		Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,50,000 /-	
32	Presto 4 microRIE (1 Unit + 1 Standard charger)	MM / LP / MP / HP / UP		Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	85,000 /-	
33	Presto 4 HP BTE (2 Units + 1 Desktop charger)			Noahlink wireless	12	Li-ion	80	30 - 110	1,25,000 /-	
34	Presto 4 HP BTE (1 Unit + 1 Desktop charger)			Noahlink wireless	12	Li-ion	80	30 - 110	70,500 /-	
35	Presto 4 NP BTE (2 Units + 1 Desktop charger)			Noahlink wireless	12	Li-ion	67	20 - 100	1,25,000 /-	
36	Presto 4 NP BTE (1 Unit + 1 Desktop charger)			Noahlink wireless	12	Li-ion	67	20 - 100	70,500 /-	
37	Presto 4 ITC/ITE (2 Units + 1 Custom charger)	LP / MP / HP / UP		Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	1,40,000 /-	
38	Presto 4 ITC/ITE (1 Unit + 1 Custom charger)	LP / MP / HP / UP		Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	80,000 /-	
39	Presto 4 ITC	LP / MP / HP / UP		Noahlink wireless	12	312	51 / 60 / 69 / 77	20 - 110	65,000 /-	
40	Presto 4 NWL CIC	LP / MP / HP		Flex cable / HI PRO	12	10A	51 / 62 / 69	20 - 105	65,000 /-	

Presto Series 4 2 Years Total Warranty

Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP in ₹
NEW		Presto CROS		Made for iPhone iPad iPod		Works with android	Made for Android™ ASHA
Software - Interton fitting 2.0 & higher		Programming Adapters - Noahlink wireless			Chipset - C 6		
41	Presto 9 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP	17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	4,05,000 /-
42	Presto 7 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP	14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,55,000 /-
43	Presto 5 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,80,000 /-
44	Presto 4 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,30,000 /-
45	Presto CROS microRIE (1 CROS" Transmitter)	LP / MP / HP		Li-ion			55,000 /-
46	Presto CROS NP BTE Transmitter (End use will pair it with a Presto 88 or 77 BTE monaural kit)			Li-ion			55,000 /-

Presto CROS 1 Year Warranty on Transmitter microRIE



	Feature Name	Move 6	Move 4	Move 3	Move 2	Move 1
Sound quality	Warp compression (WDRC), number of channels	17	12	8	6	4
	Adaptive noise reduction	✓	✓	✓	✓	✗
	Adaptive wind noise reduction	✓	✓	✓	✓	✗
	Impulse noise reduction	✓	✓	✗	✗	✗
	Microphone noise reduction	✓	✓	✓	✓	✓
	Environmental gain tuner	✓	✗	✗	✗	✗
	Environmental classifier	✓	✓	✓	✓	✓
Speech understanding	Integrated directionality	✓	✗	✗	✗	✗
	Automatic beamwidth	✓	✗	✗	✗	✗
	Combined directionality	✓	✓	✗	✗	✗
	Synchronised auto-steered directionality	✓	✓	✗	✗	✗
	Selectable beamwidth	✓	✓	✓	✗	✗
	Auto-steered directionality	✓	✓	✓	✓	✗
	Speech-focused directionality	✓	✓	✓	✓	Omni / Fixed
Feedback management	Feedback manager plus	✓	✓	✓	✓	✗
	Music mode	✓	✓	✗	✗	✗
	Preset feedback manager	✓	✓	✓	✓	✓
Apple & android connectivity / rechargeability	Direct audio streaming (MFI, android™)	✓	✓	✓	✗	✗
	Li-ion rechargeable battery	✓	✓	✓	✗	✗
Onboarding	Synchronised acclimatization manager	✓	✓	✓	✗	✗
	Acclimatization manager	✓	✓	✓	✓	✗
Convenience features	Ear to ear communication (push button, volume control)	✓	✓	✓	✗	✗
	Power-on-delay	✓	✓	✓	✓	✗
	Autophone	✓	✓	✓	✓	✓
Wireless accessories & apps	TV-Streamer 2	✓	✓	✓	✓	✗
	Remote control 2	✓	✓	✓	✓	✗
	Phone clip 2	✓	✓	✓	✓	✗
	Multi Mic	✓	✓	✓	✓	✗
	Interton Sound™ app	✓	✓	✓	✓	✓
	Remote firmware update	✓	✓	✓	✓	✓
Tinnitus solutions & Catalog	Tinnitus sound generator	✓	✓	✓	✓	✗
	Datalogging	✓	✓	✓	✓	✓

Receiver-In-the-Ear (RIE) hearing aids



MicroRIE



Rechargeable RIE



Invisible-In-the-Canal (IIC)



Completely-In-the-Canal (CIC)



Completely-In-the-Canal (CIC-W)



In-the-Canal (ITC)

Custom-made hearing aids

Behind-the-Ear (BTE) Hearing aids



SP BTE



Power BTE

Sr. No.	Model name		Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP. in ₹	
Interton Move 6						Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher						Chipset - C 5 (Coyote)			
1	Move 661 DRWC RIE (2 Units + 1 Premium Charger)	LP / MP / HP / UP		Noahlink wireless	17	Li-Ion	62 / 67 / 74 / 82	20 - 110	1,20,000 /-
2	Move 661 DRWC RIE (1 Unit + 1 Desktop Charger)	LP / MP / HP / UP		Noahlink wireless	17	Li-Ion	62/ 67 / 74 / 82	20 - 110	68,500 /-
3	Move 690 DW SP BTE			Noahlink wireless	17	675	86	60 - 120	70,000 /-
4	Move 680 DWH HP BTE Metal Hook			Noahlink wireless	17	13	80	30 - 110	52,000 /-
5	Move 670 DW BTE Thin tube / Closed NP BTE			Noahlink wireless	17	13	63 / 66	20 - 100	52,000 /-
6	Move 6 DW ITC	MP / HP / UP		Noahlink wireless	17	312	59 / 69 / 78	20 - 110	52,000 /-
7	Move 6 NWL CIC	LP / MP / HP		HI PRO	17	10A	49 / 59 / 69	20 - 105	52,000 /-
Interton Move 4						Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher						Chipset - C 5 (Coyote)			
8	Move 461 DRWC RIE (2 Units + 1 Premium charger)	LP / MP / HP / UP		Noahlink wireless	12	Li-Ion	62 / 67 / 74 / 82	20 - 110	1,04,000 /-
9	Move 461 DRWC RIE (1 Unit + 1 Desktop charger)	LP / MP / HP / UP		Noahlink wireless	12	Li-Ion	62 / 67 / 74 / 82	20 - 110	58,500 /-
10	Move 490 DW SP BTE			Noahlink wireless	12	675	86	60 - 120	49,000 /-
11	Move 461 DRW RIE	LP / MP / HP / UP		Noahlink wireless	12	312	62 / 67 / 74 / 82	20 - 110	38,000 /-
12	Move 462 DRW RIE	LP / MP / HP / UP		Noahlink wireless	12	13	62 / 67 / 74 / 82	20 - 110	38,000 /-
13	Move 480 DWH HP BTE Metal Hook			Noahlink wireless	12	13	80	30 - 110	38,000 /-
14	Move 470 DW BTE Thin tube / Closed NP BTE			Noahlink wireless	12	13	63 / 66	20 - 100	38,000 /-
15	Move 4 DW ITC	MP / HP / UP		Noahlink wireless	12	312	59 / 69 / 78	20 - 110	38,000 /-
16	Move 4 NWL CIC	LP / MP / HP		HI PRO	12	10A	49 / 59 / 69	20 - 105	38,000 /-



Rechargeable
Hearing Aids



Apple & Android
ASHA Connectivity



Water resistant
Isolate nano-coating



Interton
Sound app


Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP. in ₹	
Interton Move 3						 Made for Android™ ASHA		
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)			
17	Move 361 DRWC RIE (2 Units + 1 Desktop charger)	LP / MP / HP / UP 	Noahlink wireless	8	Li-Ion	62 / 67 / 74 / 82	20 - 110	84,000 /-
18	Move 361 DRWC RIE (1 Unit + 1 Desktop charger)	LP / MP / HP / UP 	Noahlink wireless	8	Li-Ion	62 / 67 / 74 / 82	20 - 110	48,000 /-
19	Move 390 DW SP BTE		Noahlink wireless	8	675	86	60 - 120	38,000 /-
20	Move 361 DRW RIE	LP / MP / HP / UP	Noahlink wireless	8	312	62 / 67 / 74 / 82	20 - 110	28,000 /-
21	Move 362 DRW RIE	LP / MP / HP / UP	Noahlink wireless	8	13	62 / 67 / 74 / 82	20 - 110	28,000 /-
22	Move 380 DWH HP BTE Metal hook		Noahlink wireless	8	13	80	30 - 110	28,000 /-
23	Move 370 DW BTE Thin tube / closed NP BTE		Noahlink wireless	8	13	63 / 66	20 - 100	28,000 /-
24	Move 3 DW ITC	MP / HP / UP	Noahlink wireless	8	312	59 / 69 / 78	20 - 110	28,000 /-
25	Move 3 NWL CIC	LP / MP / HP	HI PRO	8	10A	49 / 59 / 69	20 - 105	28,000 /-
Interton Move 2								
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)			
26	Move 290 DW SP BTE		Noahlink wireless	6	675	86	60 - 120	28,000 /-
27	Move 261 DRW RIE	LP / MP / HP / UP	Noahlink wireless	6	312	62 / 67 / 74 / 82	20 - 110	23,900 /-
28	Move 262 DRW RIE	LP / MP / HP / UP	Noahlink wireless	6	13	62 / 67 / 74 / 82	20 - 110	23,900 /-
29	Move 280 DWH HP BTE Metal Hook		Noahlink wireless	6	13	80	30 - 110	26,000 /-
Interton Move 1								
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)			
30	Move 180 DWH HP BTE Metal Hook		Noahlink wireless	4	13	80	30 - 110	21,000 /-
31	Move 170 DW NP BTE Thin tube / Closed		Noahlink wireless	4	13	63 / 66	20 - 100	21,000 /-

2 Years Warranty Applicable to Hearing Aid Models for which the warranty is not specified.

Note : D - Directional | R - RIE | W - Wireless | C - Chargeable



Interton Features Overview




Integrated Directionality

Adaptive microphone feature used in binaural fittings. Continuous ear-to-ear exchange of data about the environment optimises microphone settings in each ear for a wide variety of everyday listening situations.

Benefit

Users can focus on the sounds that are important to them without losing awareness of what's going on around them, so they can effortlessly follow conversations and shift attention naturally in any sound environment, all without the need for any manual adjustments.




Adaptive Noise Reduction

A unique and fast-acting spectral subtraction technology that precisely identifies speech and noise and reduces unwanted noise without degrading the speech signal.

Benefit

Noise reduction will provide comfort in noisy listening situations without compromising speech understanding.



Auto-Steered Directionality

Hearing aids analyse the sound environment and switch microphone settings between omnidirectional and directional. In directional mode, the beam width can adapt to maximise listening in noise. The switching happens so smoothly that users don't notice the change.

Benefit

Users can benefit from automatic switching for better speech understanding when needed without any need to switch programs manually. The transitions between microphone modes are performed smoothly to provide a more natural sound experience.




Adaptive Wind Noise Reduction

A specific noise reduction technique that is designed to detect and reduce wind noise only. Gain for speech sounds is preserved even in windy situations and the feature does not reduce any gain if wind is not present.

Benefit

With Adaptive Wind Noise Reduction, users won't miss out on important conversations, even in windy conditions. Whether they're cycling, hiking or golfing, they can enjoy life in any outdoor environment without compromising sound quality or speech understanding.




Speech-Focused Directionality

An adaptive directional feature, which changes directional polar patterns based on the surrounding noise source for the optimal directional benefit.

Benefit

Users can understand speech in front of them better even in dynamic sound environments.




Feedback Manager Plus

A feedback cancellation feature that uses two cancellation filters to eliminate feedback without reducing amplification or compromising sound quality. Feedback Manager Plus also has the ability to eliminate feedback from sudden movements or objects placed near the ear.

Benefit

Best-in-class feedback cancellation enables users to enjoy the sound quality without experiencing feedback or distorted sound in any situation.



Environmental Gain Tuner

Environmental Gain Tuner optimises gain settings for better listening experiences in five different environments. The hearing aids analyse the surrounding environment and automatically apply preset gain settings as the user moves from one acoustic environment to another. The gain settings can be customised for each environment based on personal preference.

Benefit

Users can benefit from automatic adjustments in gain without the need to make any adjustments in volume themselves. It's the perfect feature for active people who encounter a wide variety of listening situations. The customisable gain settings inside the feature ensure that each user experiences the best settings specific to them and the situations they encounter.




Acclimatization Manager

For some users, adapting to prescribed gain can be a challenge, especially when the users are new to hearing aids. Acclimatization Manager gradually increases gain over time to eventually achieve prescribed gain. The starting point can be set as low as 60% of the desired end gains, and the amount of time over which the gain increases can be customised to each user.

Benefit

Acclimatization Manager allows new users to adapt to amplification comfortably and seamlessly. It can also be beneficial in cases where a new user is unable to return to the dispenser for follow-up for an extended period of time.



Combined Directionality

Microphone feature used in binaural fittings. Microphone settings are set at the time of hearing aid programming by hearing care professional in each ear for a wide variety of everyday listening situations.

Benefit

Users will be able to focus on speech without losing the awareness of environment, more like a nature intended hearing.

Cable & fitting guide


	 (NOAHLINK WIRELESS)	 (HI-PRO CABLE/CS44) (ADAPTER, FLEX, 3 PIN, CS53)	 (HI-PRO CABLE/CS44) (PROGRAMMING ADAPTOR 2)	 IIC Programming	 Fitting Software
Interton Presto					Interton Fitting 2.0 or Higher
Micro RIE	●				
BTE	●				
ITC / ITE	●				
CIC		●			
IIC				●	
Interton Move					Interton Fitting 2.0 or Higher
RIE Rechargeable	●				
RIE 312 Battery	●				
RIE 13 Battery	●				
Normal Power BTE (70)	●				
High Power BTE (80)	●				
Super Power BTE (90)	●				
Custom WL	●				

How to pair your hearing aids with your Android™ device



1. Turn your hearing aids off.
2. Go to google play on your phone search for interton sound app.
3. Download interton sound app & install.
4. Go the settings menu and check that bluetooth is turned on.
follow instructions on the screen.
5. Turn your hearing aids on.
6. Click to open the interton sound app follow instructions
on the screens to connect.

How to pair your hearing aids with your iOS device



1. Turn your hearing aids off.
2. Turn on bluetooth by going to the settings menu on your iOS device.
3. Go back to the main settings menu and then, tap on general.
4. In the general menu, tap on hearing aid.
5. Turn your hearing aids on.
6. Turn your hearing aids on your iphone / ipad / ipod will detect your hearing aids.
7. Tap on pair.
8. Go to the app store, install interton sound app, follow instructions to connect to the app.

Charger feature	Premium charger (C-1)	Standard charger (C-2)	Custom charger (C-4)	Desktop charger (C-3)
Onboard power bank	YES	NO	NO	NO
LED	YES	YES limited functionality	YES limited functionality	YES limited functionality
Lid	YES	YES	YES	NO
Removable cable	USB-C	USB-C	USB-C	USB-C
Inserts (CFF)	CFF 08 for Presto MicroRIE	CFF 08 for Presto MicroRIE	Custom	CFF70 or CFF80

Premium charger



Premium Charger:

- Wireless charging solution
- Acts as powerbank
- Once charged will last for 10 Days
- Charges hearing aids for 3 times
- Fast charging enabled
- Easy to use LED indicators for battery life indication

Standard charger



Standard Charger:

- wireless charging solution
- once charged will last for 30 hours
- fast charging enabled
- comes with a lid to carry the charger

Custom charger



Custom Charger:

- Wireless charging solution
- Once charged will last for 24 hours
- Fast charging enabled
- Comes with a lid to carry the device
- The industry's only customised charger to ensure reliable charging

Desktop charger



Desktop Charger:

- Wireless charging solution
- Once charged will last 30 hrs
- Fast charging enabled

Chargers

MRP. in ₹

1	Premium Charger (Presto MicroRIE, Move - RIE 4 & 6)	15,500 /-
2	Standard Charger (Presto MicroRIE)	10,500 /-
3	Custom Charger (Presto Rechargeable ITC)	10,500 /-
4	Desktop Charger (Presto BTE, Move - RIE 3, 4 & 6)	8,500 /-

1 Year Warranty on all types of Chargers



TV-Streamer+



Multi Mic



TV-Streamer 2



Multi-Mic+

Wireless accessories

MRP. in ₹

1 TV-Streamer+ / TV-Streamer 2

17,700 /-

2 Multi-Mic+ / Multi Mic

17,700 /-

1 Year Warranty on Wireless Accessories



Discover your perfect hearing solution today



Interton GN
.....

ELITE ORGO CHEM PVT. LTD.

Sai Sangam, Office No. - 706, 7th Floor,
Plot No. 85, CBD, Belapur, Navi Mumbai - 400614.

Tel: +91-22-27564481 / 83

E-mail : eliteorgo@gmail.com

Website : www.interton.com